


REBREATHINGBAGS


■ MN 124 Disposable Rebreathingbags Tek Kullanımlık Anestezi Balonları




LATEX FREE	
CODE	SIZE
MN 124-01	0,5 LT
MN 124-02	1 LT
MN 124-03	2 LT
MN 124-04	3 LT



■ MN 125 Reusable Rebreathingbags Çok Kullanımlık Anestezi Balonları



REUSABLE BLACK RUBBER	
CODE	SIZE
MN 125-01	0,5 LT
MN 125-02	1 LT
MN 125-03	1,5 LT
MN 125-04	2 LT
MN 125-05	2,5 LT
MN 125-06	3 LT



REUSABLE SILICONE	
CODE	SIZE
MN 125-07	0,5 LT
MN 125-08	1 LT
MN 125-09	2 LT
MN 125-10	3 LT



Reusable rebreathingbags can be sterilized (autoclavable) up to 134 °C without any deformation.
Çok kullanımlık anestezi balonları 134 °C kadar otoklavda herhangi bir bozulma olmadan steril edilebilir.