

TRIBALL

Breathing Exercises System

■ MN 147 Triball-Breathing Exercises System Üç Top -Solunum Egzersiz Cihazı

TECHNICAL PROPERTIES

- Breathing exercise spirometry for pediatric and adult patients
- 600cc / s, 900cc / s and 1200cc / s flow speed ranges
- Three-dimensional chamber and color-coded breathing exercise balls
- Air filter
- Adjustable and fixed flexible hose
- Easily removable, cleanable mouthpiece
- Compact, robust and transparent tubing
- Wide base

TEKNİK ÖZELLİKLER

- *Pediyatrik ve erişkin hastalar için akciğer egzersiz spirometresi*
- *600cc/sn, 900cc/sn ve 1200cc/sn akış hız aralıkları*
- *Üç adet ölçekli odacık ve renk kodlu nefes egzersiz topları*
- *Hava filtresi*
- *Ayarlanabilir ve sabitlenebilir esnek hortum*
- *Kolayca takılıp çıkarılabilen, temizlenebilir ağızlık*
- *Kompakt, sağlam ve şeffaf gövde*
- *Geniş taban*



% 97 CALIBRATION

With the specially developed intermediary connector part measurement is made with 97% precision. Each product is tested for precision calibration.

The device consists of three graduated cylinders for indicating flow rate.

Different coloured balls in the cylinders indicate the approximate flow passing through the device.

A tube with mouth-piece is also included.

Geliştirilen özel ara konnektör ile %97 doğru ölçüm sağlanmıştır. Her ürüne hassas kalibre testi yapılmaktadır.