

■ MN 143 - 01 Valved Tee Adapter Valfli T Konnektör

TECHNICAL PROPERTIES

The Valved Tee Adapter simplifies the delivery of aerosolized medications to ventilated patients. Used in conjunction with a small volume nebulizer (SVN), the Valved Tee Adapter offers a safe and effective way to introduce the aerosol into the circuit. The Valved Tee Adapter is a permanent part of the ventilator circuit. The original spring-loaded valve opens and closes as the nebulizer is inserted and removed.

- Clean, Ready to Use
- Latex Free
- For Single Patient Use Only
- Change with the circuit

Spring loaded valve:
Opens and closes as one inserts
and removes the nebulizer.

TECHNICAL PROPERTIES

15M / 15F



■ MN 143 - 02 Inhalation Chamber with T Valved T Valfli İnhalyasyon Tedavi Çemberi

TECHNICAL PROPERTIES

JET NEBULIZER

T PIECE 10 M / 15 F

OXYGEN TUBE

